

Vantage Lite Calibration

This + means that the device needs to be calibrated.



Over time you may find that when you press keys, you don't see what you are supposed to see; or you may notice that the icons and images on the display screen don't seem to fit the display area as well as they once did. This is normal with a touch screen and can be easily fixed.

How to calibrate the screen.

Note: You must be in your CORE, the Toolbox, Page or Activity to recalibrate your screen. You cannot do this from a menu screen.

1. On the top front of the case, press and hold the **Toolbox** key while you activate the **#2 User Defined** Key. (The key that was pre-programmed as the **SET-UP** Key).
2. You now see a message that tells you to, "**Carefully press and briefly hold stylus on the center of the target. Repeat as the target moves around the screen.**"
3. With a stylus, the eraser end of a pencil or a pen cap, press on the center of the plus sign in the middle of your screen, next press on the plus in the top left corner of your screen. Press on the plus in the bottom left corner. Press on the plus in the top right corner. Press on the plus in the bottom right corner.
4. If the calibration is successful, you will see a message at the top of the screen informing you that, "**New calibration settings have been measured. Tap the screen to register all saved data. Wait for 30 seconds to cancel saved data and keep the current setting.**"
5. **Tap the screen to save your new calibration settings.**
6. Your screen should now bring you to the toolbox or the screen you started from. **If the timer continues to count down, tap the screen again until you see the screen change.** If you do not see the note, repeat these steps.

To save the Calibration Settings:

- Go to the **Toolbox**
- Then to **Maintenance Menu**
- Go to **Hardware Diagnostics**
- Select **Touch Screen** test
- Press the **TOOLBOX** key on the top of the case **5 times**. You now will hear 6 beeps in return. This verifies that your new calibration settings have been stored, so if your device loses power or is reset, you should not have to recalibrate your device.
- Press **Cancel** to exit

Liberator contacts

Regional Consultants

Mark Street: 07747016660
Verity- Jane Hart: 07879818788
Eleri Davies: 07827324286
Jacqui Malthouse: 07825512669
Paul Whitehouse: 07879418897

Office

01733 370470 (Option 2)